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Safety document for disposable mouth-nose protection mask

Product description:

This disposable mouth-nose protection mask offers effective protection for the mouth and nose area thanks to its 3-layer filter protection structure. It effectively blocks splashes and dust, which protects the health of you and those around you. The mask is made of non-woven material, filter material and ultra-fine soft fibers and offers excellent ventilation and low breathing resistance.

Important features:

- Disposable mask: the mask is for single use only and should not be reused to ensure hygienic standards.
- Low breathing resistance: The design of the mask allows easy breathing, even when worn for long periods.
- Comfortable rubber straps: These ensure a secure and comfortable fit behind the ears.
- Ergonomic fit: The mask adapts to different face shapes and can be worn by men and women.
- Sterile production environment: Guarantees a high quality and hygiene standard of the masks during production.
- Versatile: Perfectly suited for daily use in various environments where compliance with hygiene measures is required.

Instructions for use:

1. wash your hands thoroughly with soap and water before putting on the mask.
2. place the mask on the face so that it covers the mouth and nose area.
3. fasten the elastic loops behind the ears so that the mask fits securely
4. ensure that the mask fits snugly around the face and chin to avoid gaps.
5. avoid touching the mask while wearing it.
6. remove the mask by the loops after use and dispose of it immediately.
7. wash your hands thoroughly again.

Safety instructions:

- The mask should be replaced if visibly damaged or soiled.
- Keep away from children. For children, use is only recommended under the supervision of a parent or guardian.
- Not intended for reuse. Daily replacement is recommended to ensure maximum hygiene standards.

Storage:

- Store in a clean, dry place to ensure the integrity of the masks.
- Avoid direct sunlight and moisture to maintain the longevity and effectiveness of the masks.

By following these instructions and safety measures, the maximum protective effect of the disposable face mask can be guaranteed.

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